



# King of Hearts Express® AF Auto-Trigger Recorder

## Catch Events Patients Can't

The King of Hearts Express® AF monitor has the capability to record episodes of Atrial Fibrillation (AF), Tachycardia, and Bradycardia with its unique auto-triggers. In addition, the recorder may also detect certain arrhythmias that produce an *irregularly irregular* rhythm.



This strip is an example of other significant rhythms that may be automatically detected. Although the patient does not have AF, she does have a lengthening P-R interval that results in dropped beats. Untreated, the rhythm could progress to a complete heart block.

**Instromedix®**

A Card Guard Company

6779 Mesa Ridge Road, Suite 200 • San Diego, CA 92121

Phone: 800-633-3361 • Fax: 858-777-2373 • www.instromedix.com

Although the algorithm is very sophisticated, there is no guarantee that it will detect episodes of AF, nor is there any guarantee that the transition to AF will be successfully recorded in all auto-triggered events.

Copyright © 2002, Card Guard® Group of Companies. All rights reserved. Patent Pending. Card Guard® reserves the right to change specifications at any time without notice. Card Guard, LifeStar AF Express, Instromedix, Instromedix logo and King of Hearts Express are trademarks or registered trademarks of Card Guard Group of Companies.

The device that obtained this strip was a LifeStar AF Express™ recorder, which is identical to a King of Hearts Express AF recorder except for the labeling.

IMX013 Rev A