


King of Hearts Express[®] AF Auto-Trigger Recorder

Catch Events Patients Can't

The Instromedix[®] King of Hearts Express[®] AF monitor has the capability to record episodes of Atrial Fibrillation (AF), Tachycardia, and Bradycardia with its unique auto-triggers.



The following is an episode of rapid narrow complex tachycardia (PAF) that the recorder automatically captured:



Note how, in addition to the tachycardia, the recorder captured the transition from sinus rhythm to the rapid atrial activity. This provided an important diagnostic clue to the physician.

Instromedix[®]

A Card Guard Company

6779 Mesa Ridge Road, Suite 200 • San Diego, CA 92121

Phone: 800-633-3361 • Fax: 858-777-2373 • www.instromedix.com

Although the algorithm is very sophisticated, there is no guarantee that it will detect episodes of AF, nor is there any guarantee that the transition to AF will be successfully recorded in all auto-triggered events.

Copyright © 2002, Card Guard[®] Group of Companies. All rights reserved. Patent Pending. Card Guard[®] reserves the right to change specifications at any time without notice. Card Guard, LifeStar AF Express, Instromedix, Instromedix logo and King of Hearts Express are trademarks or registered trademarks of Card Guard Group of Companies.

The device that obtained this strip was a LifeStar AF Express[™] recorder, which is identical to a King of Hearts Express AF recorder except for the labeling. IMX015 Rev A