

LIFE PDA

Wireless Vital Signs Monitoring System



LifeStar™ ECG



LifeStar™ Weight Scale



LifeStar™ BP Pro



LifeStar™ Oxy Pro

The LifeStar™ Remote Monitoring System runs on the Cingular 8525, a handheld that contains all of the tools required to safely and efficiently monitor the following:

- LifeStar™ ECG
- LifeStar™ BP Pro Blood Pressure
- LifeStar™ Weight Scale
- LifeStar™ Oxy Pro Pulse Oximetry

- Portable and simple to operate
- Comfortable and rapid testing
- Anywhere, Anytime

During or after performing a test on any LifeStar™ monitor, the results are automatically transmitted via Bluetooth technology to the LifeStar™ handheld device, where it can be viewed and saved. When convenient, the data can be transmitted to the LifeStar™ Network, where the data is stored. The LifeStar™ Network provides privacy and data protection of patient information, and has the ability to track vital signs data over time.



Item #: 903-7006-00C

Reserves right to change model and/or specifications
HIPAA Compliant: HTTPD: SSL Encryption 128bit 3DES

LIFE PDA

Wireless Vital Signs Monitoring System



LifeStar™ Weight Scale



LifeStar™ ECG



LifeStar™ Network



LifeStar™ BP Pro



LifeStar™ Oxy Pro

Technical Specifications

LifeStar™ HealthKit Application

Tri-band UMTS/HSDPA; Quad-band GSM/GPRS/EDGE for global use

Microsoft® Windows Mobile™ 5.0 with Security Feature Pack

128MB ROM / 64MB SDRAM

Samsung 400MHz Processor

Sliding QWERTY keyboard

2.8" QVGA 320x240 64K Color LCD Touch Screen

Bluetooth v2.0

Integrated WiFi (802.11b/g)

Micro SD expansion slot for greater storage and expansion

* Model subject to change without notice

LIFE 
LifeWatch Technologies

LifeWatch Technologies Inc.
O'Hare International Center II
10255 West Higgins Road Suite 120
Rosemont, IL 60018
Toll Free: 800-633-3361
Fax: 800-954-2375
www.lifewatchtechnologies.com

Distributed by LifeWatch Technologies Inc. ©2007. All rights reserved. LifeStar is a trademark of LifeWatch Corp. Cingular is a registered trademark of Cingular Wireless. Windows Mobile, MSN and are trademarks or registered trademarks of MicroSoft Corporation.