

# LIFE Scale

## Wireless Weight Monitoring System

The LifeStar™ Scale is a digital weight scale enabling precise weight measurement, which can be utilized in many disease management programs or fitness regimes. An individual's weight is displayed clearly on the screen within seconds of stepping on the scale.

- Portable and simple to operate
- Step on / Step off functionality
- Self Calibration

The LifeStar™ Scale communicates with the LifeStar™ handheld wireless device and the LifeStar™ Network, provided by LifeWatch Technologies.



Item #: FG-00076C



Item #: 903-7006-00C

After weighing in, the results are automatically transmitted via Bluetooth technology to the LifeStar™ handheld device, where it can be viewed and saved. When convenient, the data can be transmitted to the LifeStar™ Network, where physicians may view patient data. Additional remote patient monitoring services are also available. The LifeStar™ Network provides privacy and data protection of patient information, and has the ability to track vital signs data over time.

# LIFE Scale

## Wireless Weight Monitoring System



Item #: FG-00076C

Technical Specifications	Parameter	Value
	Data Transmission:	Bluetooth Class I
	Transmission protocol	SPP with AT commands per ETSI TS 100,585
	Frequency	2402 ~ 2480MHz
	Modulation	GFSK
	Interchannel gap	1MH
	Power consumption	110mA Typ
	RF	100 mW
	Receive sensitivity	-80dBm Typ
	Output level	(Class1) 14dBm Typ
	Pairing	Immediate on battery insertion, pin code 6624
	Switch-on operation	Zeroed with light touch Self calibrates in seconds
	Max. stand-by duration:	Automatic shut-off after 20 seconds
	Storage capacity	20 Measurements
	Battery Type	9V Alkaline Battery life ~1000 weight procedures with Bluetooth initialization
	Operating temperature	10°C ~ 35°C
	Storage Temperature	-20°C ~ 60°C
	Dimensions (max.)	31 cm x 31 cm x 4.8 cm
	Net Weight	1063 gr
	Product specifications subject to change	



Distributed by LifeWatch Technologies Inc. ©2007. All rights reserved. LifeStar, the LifeStar symbol and HealthKit are trademarks or registered trademarks of LifeWatch Corp.

LifeWatch Technologies Inc.  
 O'Hare International Center II  
 10255 West Higgins Road Suite 120  
 Rosemont, IL 60018  
 Toll Free: 800-633-3361  
 Fax: 800-954-2375  
[www.lifewatchtechnologies.com](http://www.lifewatchtechnologies.com)